



Overview Piezomed Instruments and Sets

Bone preparation

| | | | | | | | |
|---------------------------|-----|---------|----|-----|----|------|-------|
| | | | | | | | |
| Type: | B1 | B2R/B2L | B3 | B4 | B5 | B6** | B7*** |
| Group | ③ | ③ | ③ | ③ | ③ | ③ | ③ |
| Power (Classic) ideal [%] | 100 | 80 | 45 | 100 | 45 | 100 | 100 |
| Power (Plus) ideal [%] | 85 | 60 | 35 | 85 | 35 | 85 | 85 |
| Cooling [%] | 100 | 100 | 10 | 100 | 10 | 100 | 100 |

Special lateral sinus

| | | | | | |
|---------------------------|-----|-----|-----|----|----|
| | | | | | |
| Type: | S1 | S2 | S3 | S4 | S5 |
| Group | ③ | ③ | ② | ② | ② |
| Power (Classic) ideal [%] | 100 | 45 | 100 | 80 | 80 |
| Power (Plus) ideal [%] | 85 | 35 | 100 | 80 | 80 |
| Cooling [%] | 100 | 100 | 90 | 90 | 90 |

Paro surgical periodontology

| | | |
|---------------------------|----|-----------|
| | | |
| Type: | P1 | P2RD/P2LD |
| Group | ② | ② |
| Power (Classic) ideal [%] | 80 | 80 |
| Power (Plus) ideal [%] | 80 | 80 |
| Cooling [%] | 90 | 90 |

Endo retrograde endodontics

| | | | | |
|---------------------------|-----|-----------|-----|-----------|
| | | | | |
| Type: | R1D | R2RD/R2LD | R3D | R4RD/R4LD |
| Group | ① | ① | ① | ① |
| Power (Classic) ideal [%] | 90 | 90 | 90 | 100 |
| Power (Plus) ideal [%] | 90 | 90 | 90 | 100 |
| Cooling [%] | 50 | 50 | 50 | 50 |

EX extraction

| | | |
|---------------------------|-----|-----|
| | | |
| Type: | EX1 | EX2 |
| Group | ③ | ③ |
| Power (Classic) ideal [%] | 100 | 100 |
| Power (Plus) ideal [%] | 85 | 85 |
| Cooling [%] | 100 | 100 |

Implant/ Crestal/ Anterior (A)*/ Posterior (P)* implant site preparation/ sinus lift crestal

| | | | | | | |
|---------------------------|-----|------|---------|------|---------|---------|
| | | | | | | |
| Type: | I1 | Z25P | I2A/I2P | Z35P | I3A/I3P | I4A/I4P |
| Group | ③ | ③ | ③ | ③ | ③ | ③ |
| Power (Classic) ideal [%] | 45 | 100 | 100 | 100 | 100 | 100 |
| Power (Plus) ideal [%] | 35 | 85 | 85 | 85 | 85 | 85 |
| Cooling [%] | 100 | 35 | 100 | 35 | 100 | 100 |

* Depth marking: 4 mm tip to 1st marking ring, 2 mm between markings

** Thickness saw blade: 0.25 mm, cutting depth: 8 mm

*** Thickness saw blade: 0.30 mm, cutting depth: 10 mm



Sinus lift crestal approach

I1



Z25P



I2A/I2P



Z35P



I3A / I3P



I4A/I4P



Sinus lift lateral approach


S1



S2



S3



S4



S5



Implant site preparation

I1



Z25P



I2A/I2P




Z35P



I3A / I3P



I4A/I4P



Ridge expansion

B1




B4



B6



B7



Extraction/ Explantation

EX1




EX2



B6



B7



Bone block grafting

B1



B2R/B2L



B4



B5



B6



B7



Bone chip grafting/ Bone modelling

B3



Endodontics

R1D



R2RD/R2LD



R3D



R4RD/R4LD



B3



B6




Osteotomy close to nerves

S1



S2



Corticotomy

B1



B6



B7



Periodontal surgery

P1



P2RD/P2LD



Crown preparation

S1



S2



Quick Guide Piezomed Instruments