






# assistina

## Weekly checklist





**Check liquid level!  
Refill if necessary!**



Service oil: 500 ml – 02675910	Cleaning solution: 1000 ml – 02680200
Set 02680500: 1000 ml Cleaning solution and 500 ml Service oil	





**Check function!**



**Check o-rings!  
Look for right position  
and / or  
damage!**



**Check connection is tight!**



**After consumption of one  
500 ml – bottle of Service oil  
(approx. 7000 cycles)**

**Changing of filters!**



**Once a year**