## (W&H)

### With an extra small head





The endodontics system

### Endea NiTi contra-angle handpieces for nickel-titanium files

#### Record treatment times with an extra small head

The Endea NiTi contra-angle handpieces with their extra small heads provide even better access to the treatment area.

Root canal preparation with 360° rotation is possible using both the »crown-down« and the »single-length« techniques.

Whether using an electric, an air or a torque-controlled motor, the Endea NiTi contra-angle handpiece EB-75 (16:1) enables all commercially available files to reach the required speeds.

16:1



The extremely small head gives the best possible access to the area of treatment



Maximum precision, speed and efficiency through 360° rotation

Entran cordless endodontic handpiece

Independence with an extra small head

With the cordless Entran handpiece you have absolute freedom of movement. Five different torque levels, in combination with the »auto reverse/auto forward « modes, ensure safe treatment.

The powerful, long-life Li-ion battery enables up to 30 treatments to be carried out without having to recharge.





The new, extremely small head gives the best possible access to the treatment site



Cordless endodontics with Entran

### Endea Endo Cursor contra-angle handpieces for hand files

#### Precision and safety with illumination

Use your tried and tested hand files in mechanical root canal preparation with Endea Endo Cursor. Using the »balanced force« technique, the root canal is prepared by the files' oscillating motion (approximately 60°). This is a similar motion to manual endodontics but the treatment is faster, more cost-effective and easier.





Your advantage: you can carry on using your tried and tested hand files

### Technical data

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

	Contra-angle handpieces for nickel-titanium files				Contra-angle handpieces for hand files
	Endea NiTi		Endo NiTi		Endea Endo Cursor
CE   135°C   135°C   155°C   1	Carrie		Q III	Oll more in E	
Instruments with light: Instruments without light:	EB-75	EB-79	WD-73 M	WD-74 M	EB-62
Speed reduction ratio:	16:1	2:1	70:1	128:1	4:1
Suitable inserts:	Root canal instruments suitable for 360° rotation				for root cancal
particularly nickel-titanium files, instrument sh				shaft diameter	instruments,
	2.35 mm (type 1) according to EN ISO 1797-1:1995				diameter 3.6 – 4 mm
Connection diameter:	20 mm	20 mm	20 mm	20 mm	20 mm
Coupling system:	for motors with an ISO connection according to ISO 3964 (DIN 13.940)				

# Cordless endodontics handpiece: Entran Speed: 300 rpm at the file Torque levels: 0.5 / 1.0 / 1.5 / 2.5 / 4.0 Ncm Recommended files files for rotary root canal preparation, EN ISO 1797-1:1995: shaft diameter (type 1) 2.35 mm Battery type: Li-lon, 3.7 V You can find more information at wh.com

Order Endea, Endeo and Entran from your dental supplies retailer.

Manufacturer:

**W&H Dentalwerk Bürmoos GmbH** Ignaz-Glaser-Straße 53, Postfach 1

5111 Bürmoos, Austria t + 43 6274 6236-0 f + 43 6274 6236-55 office@wh.com wh.com Distribution:

W&H (UK) Limited

6 Stroud Wood Business Centre

Park Street, St Albans, Hertfordshire AL2 2NJ, United Kingdom

t +44 1727 874990 f +44 1727 874628

f+44 1727 874628 office.uk@wh.com wh.com





